



Allevio[®] nerve block catheter

**When nothing
else works**



Allevio® nerve block catheter helps patients with migraine, Horton's disease and trigeminusneuralgy when your current treatment is not enough!

Allevio® nerve block catheter is developed in the USA where treatment is already established and used in hospitals for patients where ordinary treatment has not fully helped.

The Product

Allevio® nerve block catheter consists of two plastic catheters inserted into each other.

When used, the internal catheter is automatically bent so that it enters the back of the nose, Sphenopalatina (SPG). The purpose is to block the SPG ganglion. A 4% local anesthesia is injected through the catheter. The Allevio® nerve block catheter is CE approved as a medical device and is simple to use.

How-to?


Take a look on our webpage where we have a short video showing how the treatment is performed.

www.alleviocatheter.com



This is how to perform the treatment on a patient

1. The patient lies down on a gurney or a bed.
2. Insert the Allevio® nerve block catheter in one nostril until it stops.
3. Pull back the catheter 5 mm (a half mark, marked on the catheter).
4. Push in the inner catheter which will automatically end up in the cavity.
5. Inject 2 ml local anesthesia through the catheter.
6. Bring out the catheter.
7. Repeat the procedure in the other nostril.
8. The patient should remain in a supine position for approximately 10 minutes.



How often should treatment be performed?

Experience from hospitals in the U.S. shows that the best effect is obtained if you start with treatment two times a week during a three-week period. After the three-week period the frequency of the treatment is adjusted depending on how the patient has responded.

Some of the patients will get an immediate effect of the treatment and for some it will take a bit longer. The principle for this treatment has existed for many years, but there has not been any good method for bringing the active substance to the right place until now.

With the Allevio® nerve block catheter it is now possible to carry out the treatment in a simple way at health-care centers and clinics.

CureMed Nordic AB is a Swedish medtech company with head office and warehouse in Gothenburg. In our range, you will find world-leading products within healthcare and medical technology, many of which can be customized in production. Our goal is to always offer innovative products from leading manufacturers in a cost-effective manner to the healthcare sector.



Summary

- Now there is a possibility to help migraine patients with severe inconvenience when regular treatment is not enough
- Very easy to use
- No side affects
- CE-approved
- Used for treatment of Migraine, Hortons disease and Trigeminusneuralgy

For more information on Allevio® nerve block catheter, please contact CureMed Nordic on email: info@curemednordic.se or phone: +46 31 47 58 00

