



TYCOON
SWISS MEDICAL

Pain Management



SPECIAL CANNULAS FOR RADIOFREQUENCY AND DEVICES FOR TREATMENT OF CHRONIC PAIN

TY RaF - Radiofrequency Cannulas






TY RaF cannulas are used in nerve thermal lesions technique with the use of continuous or pulsed radiofrequency.

- Steel cannula reinforced into the tip.
- Coating made in highly insulating material.
- Large selection of sizes and lengths.
- Available in straight and curved versions.
- Marker for the ultrasound technique.
- Blunt versions to facilitate the penetration of the spinal crests and reach the vertebral foramens.
- Active tip 5 or 10 mm.
- TY RaF-TXE version, radiofrequency cannulas Slim.



Code	Needle Size		Notes	Pcs. per box
	Gauges	Length		
RF16 0100 18-R10	18 G	100 mm	Straight needle active tip 10	10
RF16 0100 18-C10	18 G	100 mm	Curved needle active tip 10	10
RF16 0100 22-R10	22 G	100 mm	Straight needle active tip 10	10
RF16 0100 22-C10	22 G	100 mm	Curved needle active tip 10	10
RF16 0145 18-C10	18 G	145 mm	Curved needle active tip 10	10
RF16 0050 18-C10	18 G	50 mm	Curved needle active tip 10	10
RF16 0100 22-R05	22 G	100 mm	Straight needle active tip 5	10
RF16 0100 22-C05	22 G	100 mm	Curved needle active tip 5	10
RF16 0050 22-R05	22 G	50 mm	Straight needle active tip 5	10
RF16 0145 22-R05	22 G	145 mm	Straight needle active tip 5	10
RF16 0145 22-R10	22 G	145 mm	Straight needle active tip 10	10

TY RaF E - Radiofrequency Echogenic Cannulas



Tycoon echogenic tip provides a top performing echo-reflection feature: all sides appear extremely bright under ultrasound. The inner diameter is not affected by the echogenic treatment, so you are free to pass wires through the cannula. This feature completely encircles the tip of the needle so all sides appear bright under ultrasound.

- Real time image.
- Faster, more accurate needle placement.
- Absence of radiation is safer for patient.
- Available in straight and curved versions.
- Available in different sizes.

Code	Needle Size		Notes	Pcs. per box
	Gauges	Length		
RF16 0100 18-R10E	18 G	100 mm	Straight needle active tip 10	10
RF16 0100 18-C10E	18 G	100 mm	Curved needle active tip 10	10
RF16 0100 22-R10E	22 G	100 mm	Straight needle active tip 10	10
RF16 0100 22-C10E	22 G	100 mm	Curved needle active tip 10	10
RF16 0145 18-C10E	18 G	145 mm	Curved needle active tip 10	10
RF16 0050 18-C10E	18 G	50 mm	Curved needle active tip 10	10
RF16 0100 22-R05E	22 G	100 mm	Straight needle active tip 5	10
RF16 0100 22-C05E	22 G	100 mm	Curved needle active tip 5	10
RF16 0050 22-R05E	22 G	50 mm	Straight needle active tip 5	10
RF16 0145 22-R05E	22 G	145 mm	Straight needle active tip 5	10
RF16 0145 22-R10E	22 G	145 mm	Straight needle active tip 10	10
RF16 1100 22-R10E	22 G	100 mm	Straight needle BLUNT active tip 10	10
RF16 1100 22-R05E	22 G	100 mm	Straight needle BLUNT active tip 5	10



V.le San Salvatore, 10 CH - 6902 Lugano/Paradiso - SWITZERLAND
telephone: +41 91 993 0504 e-mail: info@tycoonsm.com
site: www.tycoonsm.com